

Application	Amount	Volume (μL)	Conc. (ng/ μL)	Quality	Size (bp)	Comment
Illumina Short-Read Sequencing						
<i>DNA Sequencing</i>	General remarks and requirements for DNA samples: <ul style="list-style-type: none"> - Buffer: TRIS-HCL 10 mM or lowTE - 260/280 ratio 1.8-2.0 (or according to quality column) - In case of genome sequencing: high molecular for best results <p>All samples of an order must be adjusted to a uniform concentration within the specifications. Samples that do not meet the requirements listed here cannot be processed and will be rejected!</p>					
<i>Whole Genome - Low input (incl. PCR)</i>	150 – 300ng	≥ 15	> 10	see general remarks above	-	FFPE possible with constraints in output quality
<i>Whole Genome - PCR Free</i>	> 300 ng	≥ 15	> 20	see general remarks above	-	
<i>Bacterial genome sequencing</i>	15 ng	≥ 15	1	see general remarks above		
<i>Whole Exome Sequencing</i>	> 75 ng	≥ 15	5	see general remarks above	-	FFPE possible if amount based on qPCR measurement is sufficient
<i>Gene Panel</i>	600 ng	≥ 15	40	260/280 ratio 1.7-2.2 260/230 ratio 1.2	-	-
<i>Amplicons</i>	> 75 ng	≥ 15	5	260/280 ratio > 1.5	-	less material possible on request
<i>Microbiome 16S</i>	150 ng	≥ 15	10	260/280 ratio > 1.5	-	less material possible on request
<i>Microbiome Shotgun Metagenomics</i>	75 ng	≥ 15	5	260/280 ratio > 1.5	-	less material possible on request

Application	Amount	Volume (μ L)	Conc. (ng/ μ L)	Quality	Size (bp)	Comment
Illumina Short-Read Sequencing						
<i>RNA Sequencing</i>	General remarks and requirements for totalRNA samples: <ul style="list-style-type: none"> - Buffer: nuclease free water - DNase treated and cleaned up - 260/280 ratio >2.0 - RQN>=8; RQN Δ between samples <1 - ultra low input (< 1 ng total amount) possible, but must be planned beforehand <p>All samples of an order must be adjusted to a uniform concentration within the specifications. Samples that do not meet the requirements listed here cannot be processed and will be rejected!</p>					
<i>Whole Transcriptome – standard input</i>	>375 ng	>=15	25 - 80	see general remarks above	-	-
<i>Whole Transcriptome - low input</i>	>15 ng	>=15	1 - 5	see general remarks above	-	very low input amount may lead to constraints in output quality
<i>Expression Profiling (mRNA) – standard input</i>	>375 ng	>=15	25 - 80	see general remarks above	-	-
<i>Expression Profiling (mRNA) – low input</i>	>15 ng	>=15	1 - 5	see general remarks above	-	very low input amount may lead to constraints in output quality
<i>3'Prime RNA Seq -standard input</i>	>375 ng	>=15	25 - 80	see general remarks above	-	-
<i>3'Prime RNA Seq -low input</i>	>15 ng	>=15	1 - 5	see general remarks above	-	very low input amount may lead to constraints in output quality
<i>Small RNA Profiling (miRNA, Inc-RNA etc.)</i>	50 ng	>=15	50	see general remarks above	-	Low input possible on request

Application	Amount	Volume (μ L)	Conc. (ng/ μ L)	Quality	Size (bp)	Comment
<i>Sequencing of prepared NGS libraries</i>	50 ng	≥ 15	5	No primer dimer residuals	-	Buffer: TRIS-HCL 10 mM

Single-Cell RNA-Sequencing

<i>10X Genomics Chromium</i>	1.000 – 20.000 cells /sample	50 μ l	500-1.000 cells/ μ l	Vitality >80 %	-	-
<i>BD Rhapsody</i>	1.000 – 30.000 cells /sample	620 μ l	50-500 cells/ μ l	Vitality >80 %	-	-

Spatial Transcriptomics

<i>10X Genomics Visium</i>	Tissue sections placed on Visium slides, RNA extracted from tissue section RIN >7					
----------------------------	---	--	--	--	--	--

Bionano Optical Mapping

<i>Bionano Saphyr Chip</i>	1 μ g	-	>36	260/280 = 1.8 260/230 = 2.0-2.2	Megabase range	-
----------------------------	-----------	---	-----	------------------------------------	----------------	---

Application	Amount	Volume (μL)	Conc. ($\text{ng}/\mu\text{L}$)	Quality	Size (bp)	Comment
Pacific Biosciences Long-Read Sequencing						
Revio System						
<i>Transcriptome</i>						
<i>Iso-Seq</i>	>2-3 μg	>15	>200	RIN 8-10 260/280 = 2.0 260/230 = 2.2	-	DNase digested, buffer: RNase-free water
<i>Whole Genome Sequencing</i>						
<i>HiFi Reads – Standard Protocol</i>	$\geq 8 \mu\text{g}$	>120	~ 70	260/280 = 1.8 260/230 = 2.0-2.2	>50kb	homogenous HMW DNA, RNase digested, buffer: EB, DNA size distribution: 80% $\geq 50\text{Kb}$
<i>HiFi Reads – Low Input</i>	>3 μg	>50	60	260/280 = 1.8 260/230 = 2.0-2.2	>50kb	homogenous HMW DNA, RNase digested, buffer: EB, DNA size distribution: 80% $\geq 50\text{Kb}$
<i>HiFi Reads – Ultra-Low Input</i>	>50 ng	>15	>3	260/280 = 1.8 260/230 = 2.0-2.2	>50kb	homogenous HMW DNA, RNase digested, buffer: EB DNA size distribution: 80% $\geq 50\text{Kb}$
<i>Trageted Sequencing</i>						
<i>Amplicon Sequencing</i>	2-4 μg (depends on size)	50	-	-	-	clean, target-specific, buffer: EB
<i>Metagenomics</i>						
<i>Multiplexed Microbial</i>	>2-3 μg per sample	20-100	-	260/280 = 1.8 260/230 = 2.0-2.2	>30kb	Homogenous HMW DNA, RNase digested, buffer: EB DNA size distribution: 80% $\geq 30\text{Kb}$

Application	Amount	Volume (μL)	Conc. (ng/ μL)	Quality	Size (bp)	Comment
<i>Kinnex workflows</i>						
<i>Full-length RNA</i>	2-3 μg (depends on # samples to multiplex)	50	> 200	RIN 8-10 260/280 = 2.0 260/230 = 2.2	-	DNase digested, buffer: RNase-free water Up to 12-plex barcoded cDNA
<i>16S rRNA Amplicons</i>	>500 ng - 3 μg (depends on # samples to multiplex)	50	-	-	~ 1500 bp	Purified amplicons Applicants are required to order Kinnex 16S Forward and Reverse primers. These can be sourced from any oligo vendor. For optimal results, HPLC purification is recommended. Up to 384-plex samples
<i>Single-cell RNA</i>	50 ng of 10x Chromium 3' or 5' single cell DNA.	25	-	-	500-1500 bp	Applicants must provide purified Single-cell cDNA generated using the 10x Chromium Next GEM Single Cell 3' kit v3.1/ 5' kit v2 standard throughput. cDNA in the respective buffer.

Application	Amount	Volume (μL)	Conc. ($\text{ng}/\mu\text{L}$)	Quality	Size (bp)	Comment
Sequel II/IIe Systems						
<i>Transcriptome (Iso-Seq)</i>	>2 μg	>15	>130	RIN 8-10 260/280 = 2.0 260/230 = 2.2	-	DNase digested, buffer: RNase free water
<i>Whole Genome Sequencing</i>						
<i>HiFi Reads – Standard Protocol</i>	>6 μg	>120	50	260/280 = 1.8 260/230 = 2.0-2.2	>50kb	homogenous HMW DNA, RNase digested, buffer: EB DNA size distribution: 80% \geq 50Kb
<i>HiFi Reads – Low Input</i>	>1 μg	>50	20	260/280 = 1.8 260/230 = 2.0-2.2	>50kb	homogenous HMW DNA, RNase digested, buffer: EB DNA size distribution: 80% \geq 50Kb
<i>HiFi Reads – Ultra-Low Input</i>	>30 ng	>15	2	-	>50kb	homogenous HMW DNA, RNase digested, buffer: EB DNA size distribution: 80% \geq 50Kb
<i>Continous Long Reads</i>	>6 μg	>120	50	260/280 = 1.8 260/230 = 2.0-2.2	>50kb	homogenous HMW DNA, RNase digested, buffer: EB
<i>Amplicon Sequencing</i>	500 ng – 3 μg (depends on size)	50	-	-	-	clean, target-specific, buffer: EB
<i>Multiplexed Microbial</i>	>1 μg per sample	20-100	-	260/280 = 1.8 260/230 = 2.0-2.2	>30kb	-

Application	Amount	Volume (μL)	Conc. ($\text{ng}/\mu\text{L}$)	Quality	Size (bp)	Comment
Oxford Nanopore Long-Read Sequencing						
<i>Transcriptome</i>						
<i>Direct mRNA Sequencing</i>	>500 ng polyA+ or 1.5 μg Total-RNA	>12	>40	260/280 = 2.0 260/230 = 2.0-2.2	-	DNase digested, buffer: RNase free water
<i>cDNA Sequencing</i>	>200 ng polyA+ or >1 μg Total-RNA	>10	>20	260/280 = 2.0 260/230 = 2.0-2.2	-	200 ng cDNA can be used as input DNase digested, buffer: RNase free water
<i>cDNA PCR Sequencing</i>	>20 ng polyA+ or >600 ng Total-RNA	>15	>1.3	260/280 = 2.0 260/230 = 2.0-2.2	-	DNase digested, buffer: RNase free water
<i>Whole Genome Sequencing</i>						
<i>Ligation Sequencing</i>	>2 μg	>50	>40	260/280 = 1.8 260/230 = 2.0-2.2	>30 kb	homogenous HMW DNA, RNase digested, buffer: EB
<i>Ligation Sequencing with Size-Selection</i>	>10 μg	>50	>200	260/280 = 1.8 260/230 = 2.0-2.2	>30 kb	homogenous HMW DNA, RNase digested, buffer: EB
<i>Ligation Sequencing Whole-Genome Amplification</i>	>1 ng	>10	>0.1	260/280 = 1.8 260/230 = 2.0-2.2	>30 kb	homogenous HMW DNA, RNase digested, buffer: EB
<i>Native Barcoding</i>	>1 μg per sample	>50	>20	260/280 = 1.8 260/230 = 2.0-2.2	>30 kb	homogenous HMW DNA, RNase digested, buffer: EB
<i>Ultra Long Reads</i>	>50 μg	>750	>70	260/280 = 1.8 260/230 = 2.0-2.2	>100 kb	homogenous UHMW DNA, RNase digested, buffer: EEB (please inquire about EEB buffer)
<i>Rapid - gDNA</i>	>300 ng	>10	>30	260/280 = 1.8 260/230 = 2.0-2.2	>30 kb	homogenous HMW DNA, RNase digested, buffer: EB

Application	Amount	Volume (μL)	Conc. ($\text{ng}/\mu\text{L}$)	Quality	Size (bp)	Comment
Oxford Nanopore Long-Read Sequencing						
<i>Amplicon Sequencing</i>						
<i>Amplicons by Ligation</i>	>300 ng per sample	>20	>15	260/280 = 1.8 260/230 = 2.0-2.2	-	clean, target-specific, buffer: EB
<i>Rapid - 16S</i>	>20 ng	>20	>1	260/280 = 1.8 260/230 = 2.0-2.2	-	clean, target-specific, buffer: EB

Sanger Sequencing

<i>Full Service</i>						
<i>Plasmid</i>	>500 ng	-	100 – 250	260/280 = 1.8 260/230 = 2.0-2.2	<5 kb	Buffer: 10 mM Tris/HCl or water Primer concentration 10 μM
	>800 ng	-	150 - 600	260/280 = 1.8 260/230 = 2.0-2.2	5 kb – 15 kb	Buffer: 10 mM Tris/HCl or water Primer concentration 10 μM
	>1,5 μg	-	>600	260/280 = 1.8 260/230 = 2.0-2.2	>15 kb	Buffer: 10 mM Tris/HCl or water Primer concentration 10 μM
<i>Full Service / Xpress Service</i>						
<i>PCR-product</i>	max. 10 ng	-	<2,5	260/280 = 1.8 260/230 = 2.0-2.2	<100 bp	purified PCR products Primer concentration 10 μM
	50 ng	-	< 5ng	260/280 = 1.8 260/230 = 2.0-2.2	100 bp – 1 kb	purified PCR products Primer concentration 10 μM
	400 ng	-	20 - 50	260/280 = 1.8 260/230 = 2.0-2.2	>1 kb	purified PCR products Primer concentration 10 μM

Application	Amount	Volume (μL)	Conc. ($\text{ng}/\mu\text{L}$)	Quality	Size (bp)	Comment
Sanger Sequencing						
<i>Xpert Service</i>	-	-	-	-	-	Sample requirements as for Full Service
<i>Pre-mixed Service</i>						
	<i>Plasmid</i>	300 – 600 ng	-	-	-	<10 kb 5 pmol Primer/reaction 7,5 μL total volume
		>700 ng	-	-	-	>10 kb 5 pmol Primer/reaction 7,5 μL total volume
	<i>PCR-products</i>	<100 ng	-	-	-	100 bp – 1 kb 5 pmol Primer/reaction 7,5 μL total volume
		>100 ng	-	-	-	>1 kb 5 pmol Primer/reaction 7,5 μL total volume
<i>Ready-to-load Service</i>						
	<i>Plasmid, PCR-products</i>	-	20	-	-	Purified sequencing reaction

Fragment Analysis

	<i>Ready-to-run</i>	-	20	-	-	<600 bp Fully prepared samples dissolved in formamide
	<i>Pre-pared</i>	-	10	-	-	<600 bp samples dissolved in formamide without length standard
<i>STR analysis</i>		50 ng	-	5	$260/280 = 1.8$ $260/230 = 2.0-2.2$	- Buffer: 10 mM Tris/HCl, low TE or water